**75-B0024/A; 75-B0024/AZ
75-B0024/C; 75-B0024/CZ**

Centrifuga a flusso continuo

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1 Description of the apparatus

The centrifuge consist of a sheet steel case which contains the rotating spindle and the electric motor. The speed ramp is automatically controlled by an electronic circuit.

The solvent is filled through the aluminium top funnel.

This machine is proposed in two versions:

Standard model 75-B0024/A - AZ

Ventilated model 75-B0024/C - CZ

Model 75-B0024/C and 75-B0024/CZ is practically identical to the model 75-B0024/A except for the switch-rotating type which is located on top of the machine and closed in a metal box for protection from solvent vapours.

Furthermore, to avoid the excess leakage of solvent vapour from the machine and internal ventilation system has been provided for the evacuation of flammable or toxic vapour by a stainless steel cylindrical chamber ending externally to a tube 50 mm. dia. which can be connected to flexible hoses for the external evacuation or for connection to a suitable depurator.

2 Technical specification

- Electrical specs. (see label)
- Max. speed: 11.000 r.p.m.
- Centrifugal force: 4700 G
- Beaker capacity: 75 g. extracted filler approx.
- Output cap.: 15 to 30 litre/h
- Noise: 60dBA approx.
- Overall dimensions: 500 x 370 x 815 mm.
- Weight: 55 kg approx.
- Power: 550 W

3 Operating principle

This centrifuge does not require a filter, assuring therefore the highest accuracy. The test is carried out by pouring the solvent from the top funnel into the rotating aluminium beaker. Due to the centrifugal effect, the liquid spreads on the wall and moves upwards, leaving mineral particles in the beaker whilst the liquid is discharged outside through the outlet tubing. A special electronic circuit provides for a preset controlled ramp which gradually increases the speed to the maximum.

This operation takes about 10 to 15 minutes.

This centrifuge is the ideal extraction unit to complete the test started with the hot extraction apparatus which has generated a solution consisting of filler (or ash) passing the 200 mesh cloth, solvent and bitumen.

In such a situation the filterless centrifuge would required just a few minutes to complete the test.


The centrifuge is also supplied with two sieves: a top one with a wire mesh 0.146 mm. opening (no. 100 ASTM) and the second with wire mesh 0.074 mm. opening (no. 200 ASTM).

The extraction operation can also be carried out by dissolving the asphalt sample with solvent in a suitable container first and then pouring the mix into the centrifuge funnel containing these two sieves.

This method requires of a longer time than the previous one.

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
4 Set up of the machine

- 4.1 The machine has to be placed on a flat and stable working bench.
- 4.2 **⚠ Caution: The machine is supplied without oil**
- 4.3 To fill the lubricant oil supplied separately, proceed as follows:
 - 4.3.1 Remove front panel **A1** (fig. A).
 - 4.3.2 Remove stoppers **A2** pushing forward the red ring **A3** (fig. A).
 - 4.3.3 Introduce and push the metal end of the oil bottle in one of the two couplings (fig. B).
 - 4.3.4 Squeeze the bottle until the oil spurts from pos. **B2**.
 - 4.3.5 Remove the metal end of the bottle by pushing the red ring and replace the two stoppers **A2**.
 - 4.3.6 Fit again the front panel **A1**.
 - 4.3.7 Connect the centrifuge to an air depurator by the discharge tube on the rear port of the machine.
- 4.4 **⚠ Caution: Make sure that the ground is connected**
- 4.5 Move the switch on "OFF/STOP" position.
- 4.6 Check current specifications on the label before operating.
- 4.7 Connect the machine to the mains supply.
- 4.8 The apparatus has to be electrically connected to a thermal-magnetic switch (circuit breaker switch).

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5 Operating instructions

- 5.1 Remove funnel **C1** (fig. C).
- 5.2 Place completely the aluminium beaker **D1** into the rotating spindle.
- 5.3 Replace funnel **C1**.
- 5.4 Close valve **E1**, moving the little handle completely to right (fig. E).
- 5.5 Place sieves on top of funnel (if necessary).
- 5.6 Pour the solution into the sieves or directly into the funnel.

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- 5.7 ⚠ **Caution: Do not use flammable solvents**
- 5.8 ⚠ **Caution: The maximum filler content of the aluminium beaker is about 75 g.**
- 5.9 Place under the discharge tube **E2** a suitable container or bottle to collect the solution or connect a plastic or rubber tube to the tube **E2**, if it is preferred to discharge the solution far from the machine thus avoiding that the toxic solvent vapours are diffused in the laboratory. The rubber or plastic tube diameter must be big enough to avoid into the centrifuge causing serious problem to the machine.
- 5.10 Switch on the machine and wait few seconds until the maximum speed is reached.

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- 5.11 Open valve **E1** moving gradually leftward the handle (fig. E).
- 5.12 Wait until the centrifugation is completed.
- 5.13 Pour again some quantity of clean solvent to wash the funnel and other part of the machine in contact with the solution.
- 5.14 Switch off the machine and wait until the rotation is stopped.
- 5.15 Remove funnel **C1** (fig. C) and remove the aluminium beaker **D1** (fig. D).
- 5.16 Place the aluminium beaker in a oven to remove any trace of solvent and weigh the beaker to calculate the quantity of filler. Remove then the filler using a suitable flexible spatula with care to avoid damage to the beaker.
- To make easy the removal of filler from the beaker, it is suggested to line the beaker (before extraction) with a light sheet of paper that can be easily removed after the extraction together with the filler.


6 Maintenance




- 6.1 **Remove power plug of the rear panel before disassembling any part of the machine.**
- 6.2 **LUBRIFICATION**
 Use bearing oil 150 VG 10 only.
- 6.3 **CHECK OIL LEVEL**
- 6.3.1 Check oil level every 50 working hours proceeding as follows:
- 6.3.2 Remove panel **A1** (fig. A).
- 6.3.3 Remove stoppers **A2** by pushing the red ring (fig. A).
- 6.3.4 Introduce and push the metal end off bottle **B1** in one of the two couplings (fig. B).
- 6.3.5 Squeeze the bottle until the oil spurs from pos. **B2**.
- 6.3.6 Remove the metal end of the bottle by pushing the red ring and replace the two stoppers **A2**.
- 6.3.7 Fit again the front panel **A1**.
- 6.4 **OIL CHANGE**
- 6.4.1 Quantity required oil 260 ml. approx.
- 6.4.2 Change oil every 150 working hours or every two years proceeding as follows:
- 6.4.3 Remove front panel **A1** (fig. A).
- 6.4.4 Remove stoppers **A2** pushing forward the red ring **A3** (fig. A).
- 6.4.5 Remove discharge stopper placed at the bottom of the rotating spindle and the oil.
- 6.4.6 Replace discharge stopper.
- 6.4.7 Introduce and push the metal end of the oil bottle **B1** in one of the two coupling (fig. B).
- 6.4.8 Squeeze the bottle until the oil spurs from pos. **B2**.
- 6.4.9 Remove the metal end of the bottle by pushing the red ring and replace the two stoppers **A2**.
- 6.4.10 Fit again the front panel **A1**.
- 6.5 Clean periodically the external part with a suitable product.

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